# ••••• pluriSelect

# Protocol for the use of pre-filled pluriMate® centrifugation tubes

# pluriMate® - Specification

	pluriMate - 2 ml, pre-filled	pluriMate - 15 ml, pre-filled	pluriMate - 50 ml, pre-filled	
Order No. 50 pcs.	44-91002-10	44-91015-10	44-91050-10	
Order No. 100 pcs.	44-91002-11	44-91015-11	44-91050-11	
Product Description	The pluriMate <sup>®</sup> centrifugation tubes pre-filled with PBMC Spin Medium <sup>®</sup> can be used for an optimal separation peripheral blood mononuclear cells (PBMC) from whole blood and bone marrow. The key feature of pluriMate <sup>®</sup> is the porous sponge. This barrier prevents you from time-consuming and laborious overlaying of the sample material. Anticoagulated blood or bone marrow can simply be poured directly from the blood sampling tube into the pluriMate <sup>®</sup> tube. The porous barrier prevents mixture of the sample material with the separation medium. When separation is complete, the barrier prevents recontamination of the enriched cell fraction during harvest.			
Pre-filled with	PBMC Spin Medium® (Catalog 60-000	092-10)		
Enrichment of	Peripheral Blood Mononuclear Cells (F	PBMC)		
Age of blood	< 8 hours			

## Directions for the use of the pluriMate® Tube

1. Check that recommended medium, blood sample, density gradient medium and centrifuge are all at room temperature.

### Preparation of the pluriMate<sup>®</sup> Tube

2. Centrifuge at 1000 x g for 10 sec. and discard supernatant (if there is any liquid above the barrier).

#### Add Sample Material

3. Fill in sample material on top of sponge (Fig. a). Note: To reduce platelet contamination you can add pluriSpin® PLT Depletion (Order No. 19-00002-31)

	pluriMate®	pluriMate®	pluriMate®
	2 ml	15 ml	50 ml
Sample material vol.	0.25 - 1 ml	2 - 11 ml	5 - 30 ml

#### Spin

4. Centrifuge for 15 minutes at 800 x g at room temperature with in a swing bukket rotor and the **brake on.** Using blood older than 4 hours centrifuge for 30 minutes at 1000g.

#### Collect

- 5. Remove plasma by pipetting until white cell layer (Fig. d).
- 6. Collect cells in the white layer in a fresh tube (Fig. e).


#### Wash

- 7. Fill up reaction tube with wash buffer.
- 8. Spin down cells 10 minutes with 300 x g (no or small brake) at 4°C.
- 9. Pour out supernatant, leave the reaction tube on the table for 20 sec. Wash buffer excess will run down from the tube wall and collect at the bottom.
- 10. Aspirate most of the liquid above the pellet. The liquid will look foggy, these are mostly platelets aspiration will improve washing result.
- 11. Reconstitute pellet with 1 ml of wash buffer by carefully pipetting.
- 12. Repeat steps 7 to 10.
- 13. Reconstitute pellet at your desired volume.





Fig. b - Before centrifugation



#### pluriSelect - US 8912 Bancroft View Drive, San Diego (Spring Valley), CA 91977, US Phone: +1 619-202-4297 | Fax: +1 619-828-8003 support-usa@pluriselect-usa.com www.pluriselect.com