

Ministrainer 70 μ m

Catalog No	Quantity	Packaging Units
43-10070-40	25 pcs. sterile	1 Bag of 25 pcs.
43-10070-50	50 pcs. sterile	2 Bags of 25 pcs
43-10070-60	100 pcs. sterile	4 Bags of 25 pcs
43-10070-70	500 pcs. sterile	20 Bags of 25 pcs


Specification

Strainer	gray-transparent
Lot No.	see label
Mesh material	PET precision woven monofilaments, injected into the housing for mechanical stability
Housing material	LD-PE
Autoclavable	No
Chem. Stability	RTL® Buffer, TRIzol®, Isopropanol, organic solvents
To centrifuge	Yes
Shipping	Room temperature
Storage	Dry, room temperature

Component Description

The Ministrainer is a sieving device (or cell strainer) for a wide range of laboratory applications where filtration and purification of liquids is required. It has been optimized for small sample volumes. Its unique design allows fitting on a wide variation of tubes, e.g. 1.5 ml /2.0 ml reaction tubes, 15 ml conical centrifuge tubes, FACS™ tubes, Cryo vials, 24 well plate and 48 well plate.

The Pluriselect-Ministrainers are analogous or comparable with Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters.

(Legal information: Buffer RLT® is a trademark of Qiagen TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH)



Product Type	Storage Condition	Stability	Important Notes
Ministrainer	Room temperature, dry	See label	

Warning and Limitations

The product is for research and development only, not for diagnostic or therapeutic use.

Pluriselect-USA, Inc. – North America
 8912 Bancroft View Drive
 Spring Valley, CA 91977, USA
 Phone: +1 619 202-4297

support@pluriselect-usa.com
 www.pluriselect-usa.com | www.uberstrainer.com

pluriSelect Life Science, UG –worldwide
 Deutscher Platz 5e
 D-04103 Leipzig, Germany
 Phone: +49 341 333 858-99

info@pluriselect.com
 www.pluriselect.com