

Funnel

	Order No.	Size	Sterility	Packaging	Packaging Unit
<input type="checkbox"/>	42-50000-00	1 pc.	Non-sterile	Loose	/ Pack
<input type="checkbox"/>	42-50000-01	1 pc.	Sterile	Single-packed	/ Pack
<input type="checkbox"/>	42-50000-03	25 pcs.	Sterile	Single-packed	25 pcs. / Pack
<input type="checkbox"/>	42-50000-50	50 pcs.	Non-sterile	Loose	50 pcs. / Pack

Specification


LOT: stated on the label
Housing Material: PE
Sterility: Sterile




Component Description

The sterile funnel can be employed to support the sample loading onto the pluriStrainer. This way, larger sample volumes can be placed on the sieve and the maximum application volume of liquid increases from 4ml to 24 ml.

Usage of the Funnel



Attach a pluriStrainer to a sterile 50 ml centrifuge tube. Then, place the funnel on top of the strainer.



Pour sample-bead-suspension (4-24 ml) into the funnel.

Component Storage Condition

The product is stable until the expiry date stated on the label when stored unopened as indicated below. Avoid bacterial contamination.

Product Type	Storage Condition	Stability	Important Notes
Funnel	room temperature	stated on the label	-

Warning and Limitations

The product is for research and development only, not for diagnostic or therapeutic use.