

Mini-Strainer-Aid



Specification

Order No:	43-11000-50
Color:	white
Lot No.:	See label
Material:	Thermoplastic elastomer
Quantity:	50 devices, sterile
Tray material:	Polypropylene
Autoclavable:	yes (121 °C)
Storage:	Room temperature



Component Description

The Mini-Strainer-Aid is designed to facilitate the filtration of small volumes using the Pluriselect Mini-Strainer family. The perforated elastomer aid fits airtight into the Mini-Strainer. The opening holds standard 1 - 5 ml pipetted tips. The elastomer and the specially designed surface create high friction between the pipette, aid, and Mini-Strainer. That allows easy step by step handling. The Strainer-Aid reduces the filtration time in cases where centrifugation is usually needed (small mesh sizes, high suspension density, etc.).

The strainer-aid-pipette-tip-unit can be easily discarded by the tip ejection device of the pipette.

Warning and Limitation

The product is for research and development only, not for diagnostic or therapeutic use.